

## *Application Note – XpertVal 60 system*

### *XpertVal 60 cleaning procedure*

**XpertVal 60 system** - includes one complete XpertVal 60 system:

- 1 XpertVal 60 data logger
- 1 Temperature Block/Bath
- 1 Reference Probe
- 1 Computer



**The XpertVal 60 system can be cleaned with the following products:**

- **Clearklens IPA VH1**, Surface disinfectant
- **Spor-Klenz Ready to Use**, Hard Surface Antimicrobial
- **Spor-Klenz RTU**, Surface Antimicrobial

**Handling conditions of the mentioned cleaning substances:**

**Do not spill directly over the unit but only use a tissue or cleaning textile material and gently wipe the surface!**

Lives International QA  
**Ramona Kutasi**

Date: 13 August 2018

